BiliCocoon™ Disposable

BiliCocoon™ Nest Disposable

BiliCocoon™ Bag Disposable

BiliCocoon™ Bag Disposable Small

Table of contents

1. Product Description and Application.....	2
A. Description	2
B. Intended use of the Device	2
C. Additional information	2
2. Safety information	3
Types of safety precautions	3
Abbreviations and Units	6
3. Technical Features	6
4. INSTALLATION AND USE	9
A. Preliminary Precautions	9
B. Setup of the Disposable and Positioning of the Infant.....	9
C. Session settings.....	18
D. Disposal of the BiliCocoon™ Disposable	19
5. Product Catalog Numbers	19
6. Standards and guidelines	19

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User Responsibility

These Instructions for Use describe the proper setup and use of the BiliCocoon™ Disposable – the “Disposable” or the “Product”. The Product is to be used with the BiliCocoon™ Phototherapy System¹: the combination of the two of them is referred as “Device” in this document.

The Product is to be used exclusively by a properly trained user and should not be used if it is damaged, contaminated, or if parts are missing: Instead, please dispose of it and replace it.

The user is solely responsible for the risks to the patient, clinicians, third parties or properties or for a treatment with inadequate performances due to an improper use or modifications made on the Product.

Any serious incident which occurs in relation to the device should be reported to the manufacturer and the local competent authority of the state in which the user and/or patient is established. You can contact NeoMedLight at the following address: quality@neomedlight.com.

According the US federal law this Product sale is restricted to licensed medical practitioners or clinicians or under their approval.

NeoMedLight declares that this Product complies with the Regulation (EU) 2017/745.

The CE mark was obtained in 2016.

¹ Also mentioned as “BiliCocoon™ System”

1. PRODUCT DESCRIPTION AND APPLICATION

A. Description

The BiliCocoon™ Disposable is a single-use product intended to be used exclusively with the BiliCocoon™ System, a phototherapy system for treatment of unconjugated hyperbilirubinemia in the infant.

The BiliCocoon™ System consists of a light emitting Pad connected to a Light Box which emits light in the absorption spectrum of bilirubin, thus reducing the bilirubin concentration in the body of infants. *For more information, see the BiliCocoon™ System Instructions for Use.*

The BiliCocoon™ Disposable is a single-use medical device intended to be used on one patient only, thus protecting the newborn against cross-contamination.

Three sorts of Disposables are available:

- The BiliCocoon™ Nest Disposable is to be used with the BiliCocoon™ Nest System
- The BiliCocoon™ Bag Disposable is to be used with the BiliCocoon™ Bag System
- The BiliCocoon™ Bag Disposable Small is to be used with the BiliCocoon™ Bag System for small infants

For more information, see “3. Technical Features”.

B. Intended use of the Device

The BiliCocoon™ Phototherapy System is intended for the treatment of unconjugated hyperbilirubinemia² in the population of neonates and infants under 3 months old and weighing less than 10kg. It can be used in the clinical setting or in the home.

C. Additional information

- Only for the treatment of unconjugated hyperbilirubinemia, which is responsible for the majority of cases of neonatal jaundice.
- Only with medical prescription
- Only for the duration indicated in the treatment protocol
- Target patient population: neonates and infants ³.
- Over the patient's entire body except for genitals and eyes – eye protection must be worn [for the Bilicocoon Nest System](#). [Eye protection is not required for the Bilicocoon Bag System](#).
- To be used with a single patient
- To be used in hospital or at home: adapted to different environments for use (in an incubator, on a table, in a cradle, in the arms of a parent or caregiver).
- Not to be used on damaged skin

² The Device is not designed to treat cholestatic jaundice.

³ Under 3 months old and weighing less than 10 kg – definition according to EN 60601-2-50. The term “infant” used in this document is used in compliance with this definition according to EN 60601-2-50.

2. SAFETY INFORMATION

Types of safety precautions

Warning

It informs of a danger or a hazard for the infant or the operator.

Important

It mentions the actions required to obtain the expected clinical results.

Prerequisites

The Device must only be used by clinicians who are qualified and aware of the device- and treatment-associated risks and benefits; or by laypersons who are properly trained and supervised by these clinicians

Warning

Home Use: The home user must be trained before the first use of the Device about its setup and use and become aware of the device-associated risks. The home user must follow the recommendations and instructions contained in the instructions for use – concerning the System and the Disposable – and provided by clinicians

The Device is not to be used for other medical indications or on different body parts than those recommended. Improper use may lead to an ineffective treatment (i.e. no benefits) or to potential risks to the patient's health (e.g. adverse events).

Warning

Keep out of the reach of children.

Warning

Duration of Session:

Use the system only for the time prescribed. A use that does not take into account the treatment plan may lead to an ineffective treatment or to potential risks to the patient's health.

Warning

Intensive phototherapy ($>30 \mu\text{W}\cdot\text{cm}^{-2}\cdot\text{nm}^{-1}$) may not be adapted to all infants (e.g. preterm infants with a weight ≤ 1000 g)

Contraindications

The Device must not be used:

Warning

- On patients with congenital porphyria
- On patients with family history of porphyria
- In conjunction with drugs or agents that are photosensitizers
- On patients requiring a sterile environment or equipment and/or who present with skin lesions.

A use that does not take into account these contraindications may lead to an ineffective treatment (i.e. no benefits) or to potential risks to the patient's health (e.g. adverse events).

Eye Protection for Bilicocoon Nest System

Warning

Prolonged exposure to phototherapy light can cause damage to the eyes. The use of the BiliCocoon™ Phototherapy System without suitable protection of the infant's eyes may lead to risks for the infant's health, including retinal damage. Eye protection should be used for other patients who are close to the phototherapy device.

This warning is applicable only to the Bilicocoon Nest System.

Eye protection is not required for the Bilicocoon Bag System.

Patient monitoring during the phototherapy treatment

Warning

- **Bilirubin level**: The infant's bilirubin level should be regularly measured during the phototherapy treatment in accordance with the recommendations of the physician in charge.
- **Temperature**: Phototherapy can affect the patient's body temperature. The infant's temperature must be monitored according to the recommendations of the physician who is responsible for the therapy. Because the phototherapy treatment may increase the body temperature, not monitoring the infant's temperature may lead to a risk to the patient.
- **Fluid balance**: Phototherapy can affect the patient's fluid balance. Regularly monitor the patient's fluid balance and take appropriate measures to maintain this balance during phototherapy.

Warning

Side Effects for Caregivers and People near the Phototherapy System Prolonged exposure to the phototherapy device's blue light may cause discomfort to caregivers such as eye irritation, nausea, headaches or dizziness. Caregivers, staff and others who are in close proximity to the Device could be sensitive to blue light and need to protect their eyes.

Warning

Cutaneous Reactions

Cutaneous eruptions such as erythema may occur in infants treated with phototherapy.

Warning

Skin Color Changes

Blue light may hinder clinical observations by masking skin color changes, such as cyanosis.

Warning

Rare Allergic Reactions

Despite the biocompatibility evaluation tests conducted in accordance with the state of the art, rare allergic reactions may nevertheless occur.

Warning

Cross Contamination:

The Pad is to be used in combination with its specific Disposable. The infant should not be placed in direct contact with the Pad.

The Disposable is not designed for re-use and should not be cleaned, disinfected, or placed in contact with the cleaning products used to clean the Pad.

Warning

Not using the Disposable or having the cleaning products in contact with the infant may lead to allergic reactions for the infant.

The BiliCocoon™ Disposable is not intended to be used with a phototherapy treatment system other than the BiliCocoon™.

The Disposable is a **single use, single-patient** medical device. It should be used respecting the usual hygiene precautions and only use this single-use device for one patient. There is indeed a cross-contamination risk if the Disposable is used for more than one patient.

Disposable may be used multiples times on the same over a period of 24 hours. However, the user must make sure that the Disposable is neither damaged nor soiled between each use

Warning

Variation of the Ambient Conditions

Variations of the ambient conditions (room temperature and humidity, sun exposure, nearby devices) may unfavorably affect the patient's health, including their temperature or their fluid balance.

Moreover, use of this Device outside of the given operating conditions may compromise its functionality.

Warning

Reflective Foils

Never use reflective foils to increase the effectiveness of phototherapy treatment, it may hazardously increase the body temperature.

Warning

Flammability

Keep away from flame, sparks and other ignition sources.

Warning

No modification of the Disposable is allowed.

The use of a Product which has been modified disregarding the manufacturer's indications may lead to risks to the patient (e.g. adverse events), the user, other peoples or properties.

Warning

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Do not place the infant inside the Disposable.


For all the phototherapy-related risks refer to the BiliCocoon™ System Instructions for Use.




Warning

Labels Integrity and Readability

In case labels are degraded, the information addressing the safety of the patient, users, third parties, and properties can be found in this document. If the product cannot be identified by its label because the latter is damaged, please use another item.

Symbol	Description	Symbol	Description
	Manufacturer		CE Marking declaration according to Regulation (EU) 2017/745
	Consult instructions for use		Date of manufacture
	Catalogue number		Batch code
	Protect infants' eyes with opaque eye protection (Bilicocoon Nest)		Do not re-use
	Do not use if the package is damaged		Non-sterile
	Absence of Latex		Operating temperatures limits and storage temperature limits for products out of the box

 Storage temperatures limits for Products in the original box	 Keep away from sunlight
 Keep dry	 Medical device
 Position of the infant in the Bag Disposable and in the Bag Disposable Small	 Position of the infant in the Nest Disposable

Abbreviations and Units

Abbreviations	Units and Description	
°C	Degree Celsius	(Unit of Temperature)
kg	Kilogram	(Unit of Mass)
m	Meter	(Unit of length)
hPa	Hectopascal	(Unit of atmospheric pressure)


3. TECHNICAL FEATURES

The BiliCocoon™ Disposable is a non-woven medical device intended to be used exclusively with the BiliCocoon™ Pad, in order to keep the phototherapy system clean during treatment (*see the BiliCocoon™ System Instructions for Use*) and provide a comfortable interface between the skin of the baby and the Pad.

The BiliCocoon™ Disposable is a single-use product, intended to be used for one patient only and to prevent cross contaminations.


- The BiliCocoon™ Nest Disposable is to be used with the BiliCocoon™ Nest System
- The BiliCocoon™ Bag Disposable is to be used with the BiliCocoon™ Bag System
- The BiliCocoon™ Bag Disposable Small is to be used with the BiliCocoon™ Bag System for small infants

BiliCocoon™ Nest Pad

1. Lighting Surface
2. Cable
3. Connection

BiliCocoon™ Bag Pad

1. Lighting Surfaces
2. Cable
3. Connection

- **Nest Disposable**

Front view
Contact with the patient


Insertion slit of the
Nest Pad

Back view

- **Bag Disposable / Bag Disposable Small**

The Bag Disposable and the Bag Disposable Small are similar regarding their design. They only differ in the position of the infant fastening system, which is slightly higher in the Bag Disposable Small.

Use & Storage (For Products out of the box) Conditions

Temperature	10°C – 35°C ⁴
Relative humidity	15% – 90% without condensation
Atmospheric pressure	700 – 1060 hPa ⁵
Others	Indoor use & storage only

Disregarding the operating and storage conditions may lead to the Product deterioration, therefore creating risks for the patient, the users, third parties or properties

⁴ Equivalent to 50-95°F

⁵ Equivalent altitude: 0 – 3000 m / 0 – 9842,52ft

4. INSTALLATION AND USE

A. Preliminary Precautions

Select the Disposable according to the Pad type (Nest Pad or Bag Pad).

Important

If a Bag Pad is used, a Bag Disposable or a Bag Disposable Small should be used depending on the size of the infant. For smaller babies, it is recommended to use a Bag Disposable Small.

Warning

Before inserting the Pad in the Disposable, check the Pad mean irradiance as described in the protocol contained in the BiliCocoon™ System Instructions for Use.

If the device is not verified before each use, there is a risk of under-exposure – leading to an ineffective treatment – or of over-exposure.

Warning

Make sure that the enclosure of the BiliCocoon™ Pad is neither dirty nor broken, and that there are no sharp edges that might cause an injury. Clean the Pad before each use as described in the BiliCocoon™ System Instructions for Use.

Warning

Ensure that the Pad is clean and dry. The BiliCocoon™ Disposable should not come in contact with the detergents used to clean the Pad.

Warning

Never place the infant directly on the Pad. Always use the Disposable specifically designed for the BiliCocoon™ System.

Warning

Infants should wear a diaper during the phototherapy session.

Warning

Check that the packaging of the BiliCocoon™ Disposable is not damaged: The use of a Product whose cleanliness is intended as a main feature and controlled thereof, but whose packaging is damaged, may lead to a risk for the patient's health.

Warning

Check that the Disposable is clean and not damaged. The use of a Disposable that is damaged or contaminated, use of the incorrect disposable, or the use of the BiliCocoon™ System without Disposable may lead to an inadequate treatment or to risk for the patient.

Warning

Ensure that the fastening tabs, their adhesion surfaces, the Pad fixing tabs, and the Infant fastening system are intact.

Warning

The phototherapy system may increase the patient's body temperature when it is used in combination with a thermotherapy system (e.g. infant incubators, infant transport incubators, infant radiant warmers, devices supplying heat via blankets, pads or mattresses). In this case, the user should measure the infant's body temperature directly (i.e. on the skin's surface) and use the baby controlled mode of these devices, otherwise the set air temperature of the incubator or the heater output of the radiant warmer or heated mattress has to be adjusted according to the body temperature measurements.

Place the Light Box outside of these systems.

B. Setup of the Disposable and Positioning of the Infant

Warning


Never use the BiliCocoon™ System without the associated Disposable.

Important

We recommend not to use the Bag System for babies weighing less than 2500g


- For the BiliCocoon™ Nest Disposable

Front View


1. Soft part
2. Strap used to keep the Pad bent over the newborn
3. Insertion slit of the Pad

Back View

4. Back face
 5. Side tab
 6. Adhesion band
- The side tab (5) is meant to be fixed to the adhesion band (6)*

1. Place the BiliCocoon™ Disposable on a flat surface with the soft part (1) facing the user.

2. Insert the Pad in the Disposable through the insertion slit (3) so that:

- the Disposable completely covers the Pad
- The Pad Lighting Surface (the one with the label) is immediately under the aforementioned soft part (1)

Bend the Pad along the direction of the optical fibers (as in the picture) to insert it more easily.

Warning 3. Check that this soft part (1) is in contact with the Pad Lighting Surface. The Lighting Surface is not to be placed in contact with the Disposable's back face (4).

Important Disposable.

4. Check that the Pad is fully inserted to the edge of the BiliCocoon™ Disposable.

Warning 5. After the Pad has been correctly inserted into the BiliCocoon™ Disposable, check that the Disposable has not been damaged during the insertion. A damaged Disposable may lead to an inadequate treatment or a risk to the patient

Warning 6. The infant must wear a diaper during treatment. Place the diaper trying to maximize the skin surface area which is in contact with the Lighting Surface.

7. Place the infant on the Disposable with the baby's legs towards the Pad cable.


Warning

Never place the infant's head close to the Pad cable. This may otherwise lead to an ineffective treatment or to a risk of strangulation for the patient

Warning 8. Protect infant eyes with suitable eye protection for BiliCocoon Nest System. Eye protection is not mandatory for the BiliCocoon Bag System.


9. The user may keep the *BiliCocoon™ Nest Pad* sides bent over the infant through the appropriate strap (2). Fix the side tab (5) on the adhesion band (6). This disposition let the blue light come from the sides with a larger body surface which is exposed.

10. Start the treatment


- For the BiliCocoon™ Bag Disposable and the BiliCocoon™ Bag Disposable Small

**Folded Bag Disposable/
Bag Disposable Small**


1. Front part with the patterned adhesion surface
2. Insertion slits (x2) for the two Pad Lighting Surfaces
3. Pad fixing tabs

Unfolded Bag Disposable/Bag Disposable Small

4. Lateral flaps with fastening tabs
5. Shoulder fastening tabs
6. Infant fastening system

Infant Fastening system


7. Side tabs
8. Ventral flap

1. Place the BiliCocoon™ Disposable on a flat surface

2. Place the Disposable on a table as in the picture, i.e. with the front part (1) facing the user when the Disposable is folded.


4. Insert the second Pad Lighting Surface in the front part insertion slit (indicated by the arrow in the picture).


The Pad's Lighting Surfaces should be completely covered by the Disposable and should be in contact with the soft inner surface of the BiliCocoon™ Disposable.

3. Insert the Pad Lighting Surface with the label in the insertion slit below (indicated by the arrow in the picture).


Warning: Never place the Pad Lighting Surface with the label in the insertion slit of the disposable's front part. Disregarding this instruction may lead to a patient with the respiratory tract constricted.

Important 5. Check that the two Pad's Lighting Surfaces are fully inserted to the edge of the BiliCocoon™ Disposable.


6. Cross and secure the Pad fixing tabs (3) so that the Pad cannot slip out of the Disposable


Warning It is important to cross these tabs as shown in the picture on the right **prior to placing the infant in the Disposable**.

7. Unfold the disposable and its infant fastening system (6).

Warning 8. After the Pad has been correctly inserted in the Disposable, check that the Disposable has not been damaged during the insertion. A damaged Disposable may lead to an inadequate treatment or to a risk to the patient.

9. Place the infant on the Disposable, using the neckline cut-out as a reference to adjust his/her position.

Warning The infant must wear a diaper. Place the diaper trying to maximize the skin surface area which is in contact with the Lighting Surface.

10. Check that the infant is in contact with the Disposable inner surfaces (soft face).

11. Place the ventral flap (8) between the infant's legs and fold it on his/her stomach. Attach the side tabs (7) on the ventral flap (8).

Important 12. If the ventral flap (8) is too long, fold it in order to have the stomach exposed and to maximize the body surface area that will be treated.


Important 13. Check that the infant is properly positioned and securely fastened in the Disposable.

14. Fold the front part (1) on the infant chest so that the patterned adhesion surface faces the user.

Warning: Make sure that the infant's head is outside of the Bag.


15. Position the shoulder fastening tabs (5) on the patterned adhesion surface (1). This will prevent light emission outside of the Disposable.

Warning: Be careful not to bother the respiratory tract. It is important not to cross the fastening tabs.

16. Fold the lateral flaps (4) to wrap the infant and prevent light emission outside of the Disposable. Place the fastening tabs of the lateral flaps (4) on the patterned adhesion surface (1).

17. Eye protection is not required for the Bilicocoon Bag System.

18. Start the treatment.

NB: For a more comfortable treatment, it is possible to open the upper corners and let the infant's arms pass through them.

NB: A Disposable can be reused – on the same patient exclusively – after the end of the session and over a period of 24 hours only if it's clean and undamaged. Between uses, fold the Velcro adhesive tabs on themselves to preserve their adhesiveness.

C. Session settings

Refer to the BiliCocoon™ System Instructions for Use.

D. Disposal of the BiliCocoon™ Disposable

Warning

The BiliCocoon™ Disposable may be soiled during a phototherapy session. It must be disposed of according to local practices for contaminated products.

Warning

Contaminated products are to be disposed of separately from consumer waste. Contaminated products may otherwise harm people and pollute the environment.

The BiliCocoon™ Disposable is not to be gas-sterilized or autoclaved.

5. PRODUCT CATALOG NUMBERS

BCC010	BiliCocoon™ Light Box
BCC101	BiliCocoon™ Nest System
BCC111	BiliCocoon™ Nest Pad
BCC12150	BiliCocoon™ Nest Disposable x50
BCC201	BiliCocoon™ Bag System
BCC211	BiliCocoon™ Bag Pad
BCC22150	BiliCocoon™ Bag Disposable x50
BCC22150-S	BiliCocoon™ Bag Disposable Small x50 (Not available in the United States)
BCC301	BiliCocoon™ Fixation System

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